

**DISRUPT
CAD III**

**MAKING
BIG WAVES
IN COMPLEX CALCIUM**

**Largest and
Most Rigorous
Shockwave Study**

SHOWCASES SAFETY, EFFECTIVENESS & EASE OF USE OF CORONARY IVL

CAD III By The Numbers

384

Patients at
47 Sites

100%

Severe Ca++

47.9mm

Calcium Length

98%

IVL Crossing &
Therapy Delivery

99%

Stent Delivery

0.3%

Final Major
Dissections

0.3%

Final Perforations

0.3%

Abrupt Closure

0%

Slow Flow/
No Reflow

11.9%

Residual
Stenosis

1.7 mm

Acute Gain

Presented at TCT Connect and Simultaneously Published in JACC: J. Hill, D. Kereiakes, et al. "Intravascular Lithotripsy for Treatment of Severely Calcified Coronary Artery Disease: The Disrupt CAD III Study." JACC 2020.

Key Findings

CAD III confirms IVL safety with low rates of major peri-procedural clinical and angiographic complications, setting a new bar for safety in complex coronary lesions.

- 92.2% of patients were free from MACE at 30 days, a composite of CD (0.5%), MI (7.3%), or TVR (1.6%)
 - Primary MACE driver was in-hospital non-Q-wave MI (5.7%)
- Low risk of perforation (0.3%), major dissection (0.3%), abrupt closure (0.3%), and slow flow/no reflow (0.0%) at the end of procedure

CAD III showcases IVL effectiveness with large lumen gains that facilitate stent delivery and optimize stent expansion.

- 92.4% procedural success rate, defined as successful stent delivery (99%) residual stenosis <50% (100%), & without in-hospital MACE (93%)
- Successful IVL crossing & therapy delivery in 98% of lesions, correlating to 99% stent delivery
- 1.7mm acute gain and 11.9% final in-stent residual stenosis

CAD III demonstrated Coronary IVL's ease of use and quick learning curve to achieve consistently predictable outcomes.

- Despite >80% of operators having no prior experience with IVL, MACE, procedural success and device crossing success were similar between first roll-in first case and pivotal cohort



DISRUPT CAD III Design & Patient Demographics

Stable angina, unstable angina or silent ischemia



Heavily calcified, de novo coronary lesions
RVD 2.5 – 4.0 mm, stenosis ≥50%,
Lesion length ≤40 mm

Objective: Prospective, multicenter, single-arm global IDE to evaluate the safety and effectiveness of coronary IVL

Primary Safety Endpoint: Freedom from MACE at 30 days

- Cardiac death, Myocardial infarction, or Target vessel revascularization

Primary Effectiveness Endpoint: Procedural success


- Successful stent delivery with residual stenosis <50% and without in-hospital MACE

Secondary Performance Endpoints: Clinical and Angiographic Success


Baseline Characteristics:

71 yo	65% DS	57% LAD, 13% LCX, 29% RCA, 1% LM
40% DM	48mm Ca++ length	292 Ca++Arc @ Max Ca++ Site
3.0mm RVD	25mm Lesion Length	.96mm Thick Ca++ @ Max Ca++ Site

Large Circumferential Lumen Gains to Facilitate Stent Delivery, Apposition & Expansion



Low Rates of Complications, Similar to DISRUPT CAD I & II



Low Rates of MACE, Mostly Driven by In-Hospital NQWMI

